

PRODUCT DESCRIPTION

- INSTOP is a brand name for Heparin Cap.
- It is also referred as “injection stopper” or as “yellow stopper” used to close the female luer of any compatible device when not in use but still in place at the patient end.
- Designed for providing access to the catheter for injecting fluids or drawing blood through the injection membrane, using syringe and the hypodermic needle.
- Injection membrane tested as per ISO 8536-4 for self-sealing.
- Luer conforming to ISO 80369-7 for suiting the female end of any standard product.
- INSTOP are made of a variety of high-quality materials, such high-quality medical-grade silicone rubber used to administer sterile injectable fluids such as pharmaceutical agents.
- Easy to install and uninstall. Able to connect with a wide range of connectors.
- Prevents cross-contamination infections and increases security during infusions.
- INSTOP is sterile-packed individually using Ethylene Oxide (EO) gas according to ISO 11135.
- The products are packaged with blister packs with a definitive usage life-span marked on the packs. Blister packs provide protective barrier that helps safeguard products from external factors.
- Registered with Medical Device Authority of Malaysia under Act 737.



DATASHEET

PRODUCT CODE	PACKAGING				
	PACK (PCS)	CARTON (PACK)	PACK DIMENSION L x W x H (CM)	CARTON DIMENSION L x W x H (CM)	PACK WEIGHT (KG)
100U - 11	100	50	11.0 x 11.0 x 6.0	56.5 x 23.5 x 31.0	0.15 +-

Distributed by:

Intaniaga Sdn. Bhd. (658891-X)

Lot 2, Jalan CJ 1/4, Kawasan Perusahaan Cheras Jaya, 43200 Cheras, Selangor, Malaysia.

+603-9076 2010 / 3020

sales_isb@muzamal.my

www.muzamal.my

MDA Reg. No.:
GMD22056769318A

Advertisement No.:
MDAMD 0027/2024

Ref. Doc. No.:
ISB-ART-SMD-07 (REV00)

Effective Date:
1 MAY 2024

